SE Manual för Reläbox

GB User Manual Relay Unit GL-RJ

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
GB User Manual Relay Unit GL-RJ

General

GL-RJ is a relay unit used to wirelessly switch on and off lamps, tape recorders, radio, TV sets, computers and electric heaters. The relay unit is equipped with built-in IR-receiver for control by a GewaLink transmitter. An overload protection is built-in.

The relay unit can be programmed to work in four different ways:


- To switch off and on by the same button on the transmitter.
- To switch on by one button and switch off by another button on the transmitter.
- To switch on and switch off by the same button on the transmitter but also a specific button to switch off. To be used when several devices in the same room are to be switched off by one and the same button.
- Monostable function, which means that the outlets on the relay unit are power supplied as soon as the button on the transmitter is pressed.

Connection

1. Connect the relay unit to a wall outlet.
2. Connect the equipment, which is to be controlled, to the outlet on the relay unit.
3. Check that the equipment is power supplied by pressing the programming buttons for the two left and the two right outlets.

Programming buttons for the two left and the two right outlets. Can also be used to control the outlets without an IR transmitter.

Programming

The relay unit can be programmed by the GewaLink channels 00-61, 63-127 and coded channels (4096). There are four different ways to program the relay unit. The first alternative is most frequently used.

Alt. 1 Programming "On/Off" by the same button on the transmitter

Press the programming button representing the outlets you want to control on the relay unit and press, at the same time, any button on the IR transmitter until the relay clicks 2 (3) times.

Alt. 2 Programming "On/Off" with two different buttons on the transmitter

Press the programming button representing the outlets you want to control on the relay unit and press, at the same time, any button on the IR transmitter until the relay clicks 2 (3) times, change to another button ("Off") on the IR transmitter until the relay clicks 2+2 times.

Alt. 3 Programming "On/Off" with one and the same button and another button for "Off" on the transmitter

This method makes it possible to switch on/off each outlet and switch off all outlets in the same room simultaneously. By this procedure "Off" has to be programmed on all outlets with the same channel.

Press the programming button for the outlets you want to control on the relay unit and press at the same time any button on the IR transmitter until the relay clicks 2 (3) times. Change to another button ("Off") on the IR transmitter until the relay clicks 2 times. Then program the next outlet with another channel but with the same channel used for "Off".

Alt. 4 Programming Monostable Function

By this procedure power is supplied as long as the button on the transmitter is pressed. When the button is released the current supply stops.

Press the programming button for the outlets you want to control on the relay unit and press at the same time any button the IR transmitter until the relay clicks 2 (3) + 2 times.

Overload protection

If GL-RJ is overloaded the overload protection will be activated. Press the button for the overload protection to reset, see figure.

Cleaning

When being cleaned the relay unit must not be connected to the mains. To clean the relay unit, wet a cloth in a weak solution of washing up liquid and water, wring out thoroughly and wipe.

Technical data

Voltage:	230V AC, 50 Hz
Max load:	8 A
Power consumption (stand-by):	Ca 10 mA stand-by
Programming channels:	GewaLink channels 00-61, 63-127 and 4096-code
Surroundings:	Indoors

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Gewa AB