- SE** **Manual för Access USB/Keyboard Interface**
- GB** **Manual for Access USB/Keyboard Interface**
- FR** **Manuel d'utilisation en Français
Interface Access USB Clavier /
Souris**
- DE** **Gebrauchsanleitung für Access
USB/Keyboard Interface**

Innehåll, Contents, Table des matières, Inhalt

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GB Manual for Access USB/Keyboard Interface

Introduction

ACCESS USB/KEYBOARD INTERFACE enables wireless data communications with infrared light to a personal computer. The interface is a receiver that takes in IR signals and converts them so that the computer can understand what has been transmitted. Data communications can be for both keyboard and mouse functions.

Interface Description

Transfer of keyboard data is accomplished with Gewa's special IR codes for keyboard and mouse data. These codes can be generated by the following IR transmitter models:

Gewa Prog III Module, Gewa Prog III Micro Module, and the Access Medi and Access Maxi keyboards.

There are normally two ways of controlling mouse functions on a computer:

1. By moving a conventional mouse in different directions.
2. By pressing the keys on the numeric section of the computer's keyboard.
(NOTE: Mouse keys must be activated by Windows settings.)

Depending on how the interface is connected to the computer, it can simulate mouse control in one (via a keyboard port) or both (via a USB port) of these ways.

	<p>Simulating conventional mouse control</p>  <p>The following transmitters can simulate conventional mouse control: GP III Module GP III Micro Module</p>	<p>Simulating numeric mouse control</p>  <p>The following transmitters can simulate numeric mouse control: GP III Module GP III Micro Module Access Medi keyboard Access Maxi keyboard</p>	<p>Simulating other keyboard data:</p>  <p>The following transmitters can simulate keyboard data: GP III Module GP III Micro Module Access Medi keyboard Access Maxi keyboard</p>
Computer's USB port	Yes	Yes	Yes
Computer's keyboard port	No	Yes	Yes

Connection

If the computer has full support for USB, it is always best to connect the interface to the USB port.

Connection of Interface via a USB Port

For the interface to function as a USB interface, the computer must have full support for USB, which is integrated into Windows 98 and later Windows versions. USB support can be added to previous versions, but new software may need to be installed. Contact a Microsoft retailer for more information.


The interface may also be connected to Macintosh computers that have support for USB.

When the USB port is used, the computer permits several keyboards and mice to be used at the same time.

To connect the interface:

1. Make sure that the computer is turned off.

2. Connect the provided USB cable between the interface and the computer's USB port. No other cables shall be connected.
3. Start the computer. A dialog box may appear prompting you to insert the operating system's installation CD.
4. Select the same ID as for the transmitter unit. See Installing the ID Selector for more information.
5. Place the interface at the computer, making sure that the IR lens at the front of the unit is unobstructed.


Connecting the Interface via a Keyboard Port

For a computer without USB ports, the interface is connected to the computer's keyboard port, which means that the computer's standard keyboard cannot be connected as usual. To solve this problem, a branching adapter can be used. The branching adapter is not provided and must be purchased separately.

To connect the interface:

1. Make sure that the computer is turned off.
2. Connect the provided keyboard cable between the interface's keyboard port and the computer's keyboard port. If the cable connector doesn't fit the computer's keyboard port, an adapter must be connected between the connector and the computer. This adapter is not provided.
If a conventional keyboard will be used at the same time, a branching adapter must be used as shown in the figure.
3. Start the computer.
4. Set the same ID as for the transmitter unit. See Installing the ID Selector for more information.
5. Place the interface at the computer, making sure that the IR lens at the front of the unit is unobstructed.

Installing the ID Selector


The ID selector is used to identify the receiver and transmitter units. The ID selector has 16 different positions, which means that several transmitter units can be operated in the same room without interfering with one another. If the receiver's ID is set to position 0, all of the various IDs can be received simultaneously. The other positions require that the transmitter ID be in agreement with the receiver ID.

Indicator Lamps

The interface is equipped with four red indicator lamps mounted behind a red plastic window. The lamps may be on, off or blink, depending on current activity.

- Num Lock, Caps Lock and Scroll Lock: On when these functions are active.
- Data: Blinks when keyboard or mouse data is received by the interface.
- All lamps blink in unison:
 1. At power-up.
 2. If a faulty ID code is selected.
 3. For unknown data.

Interface from the front

Miscellaneous

Cleaning

Make sure that the computer is turned off. Clean the interface with **very hard-wrung cloth** that has been dipped in water with a small amount of dishwashing detergent. Wipe with a dry cloth after cleaning.

Technical Data

Internal power supply:	From the computer's USB port or keyboard port
Ports:	USB and keyboard (PS-2 standard)
Transmission, mouse data:	Gewa mouse standard
Transmission, keyboard data:	Gewa keyboard standard
EMC standard emission:	EN 50081-1
EMC standard immunity:	EN 50082-1
Temperature range:	0-40°C
Size, W x D x H	95 x 55 x 23 mm
Humidity and dust protection:	IP-20

Included Accessories

Art. no.	Designation	Description	Total
663210	USB cable	USB connectors, 1.5 m	1
663214	Keyboard cable	6-pole mini DIN, 1.5 m	1

Optional Accessories

Art. no.	Designation	Descriptions
6556	Keyboard adapter	6-pole mini DIN/5-pole DIN
663212	Branching adapter, keyboard	Branching adapter with 6-pole mini DIN

Compliance with Medical Device Directive

Access USB/Keyboard Interface is marketed as a technical aid for persons with disabilities. It complies with the Medical Device Directive 93/42/EEC, which includes, among other things, applicable requirements for EMC and electrical safety.

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Wir als Hersteller: GEWA AB
Gewa AB